

Purpose

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The PAI-1 ELISA Controls have been developed to provide customers with control material for quality monitoring of day-to-day assay performance. The expected recovery ranges for each of these controls are listed below.

Disclaimer

Specific claims as to the expected ranges of PAI-1 levels found in plasma are not made by the manufacturer, and must be determined by independent survey of human plasma specimens.

Contents

Each package contains the following materials.

Precision Controls (one ea): PAI-1 in Buffer, at 1.0 mL.

Vial	Part number	Lot number	Dosage
Level A	17-7022A	7040	High
Level B	17-7022B	7040	Mid
Level C	17-7022C	7040	Low

Use

The PAI-1 ELISA Controls are supplied ready-to-use. Controls should be added to the PAI-1 ELISA microtiter plate at 100 µl/well and the ELISA procedure followed as detailed in the Directions for Use booklet which accompanies each kit.

Storage

Store undiluted PAI-1 ELISA Controls at 2° - 8°C. Protect from light.

Data PAI-1 ELISA Control Reference Ranges:

Control #	Mean ng/mL	ng/mL Range
Control A	1.15	1.00 – 1.30
Control B	0.80	0.70 – 0.90
Control C	0.19	0.17 – 0.22

Manufactured By:

Siemens Medical Solutions Diagnostics
Tarrytown, NY 10591-5097 USA

Technical Support:

Siemens Medical Solutions Diagnostics
Technical Care Center

In the United States:

tel: (877)229-3711

For all other areas:

tel: +1 781 551 7591

email: tcconco gene.med@siemens.com