HER-2/neu ELISA Controls REF 06489884

INSTRUCTIONS FOR USE

Purpose

FOR IN VITRO DIAGNOSTIC USE

The HER-2/neu ELISA Controls have been developed to provide customers with control material for quality monitoring of day-to-day assay performance. The expected recovery ranges for each of these controls are listed below.

Contents

Each package contains the following materials.

HER-2/neu p105 ECD in buffer with 0.09% sodium azide, vialèd at 0.5 mL each.

This product is prepared from material produced by recombinant methods with non-human cell lines.

Vial	Part number	Lot number	Level
A	10-7022A	7043	High
B	10-7022B	7043	Mid
C	10-7022C	7043	Low

Use

The HER-2/neu ELISA Controls are treated just as a sample prior to use in the HER-2/neu ELISA. For sample preparation methods, and for technical assistance, refer to the ELISA IFU, REF 06489876.

CAUTION: This device contains material of human or animal origin and should be handled as a potential carrier and transmitter of disease.

Storage

Store undiluted HER-2/neu ELISA Controls at 2–8°C. Do not use beyond the expiration date indicated. Protect from light.

Controls that have been diluted in Sample Diluent may be stored at 2–8°C, protected from light, and used for up to 2 weeks.

Data HER-2/neu ELISA Control Reference Ranges. Expected values are traceable to the standardization of the HER-2/neu ELISA. Refer to the ELISA IFU for further details (REF 06489876).

Control #	Dilution	Mean ng/mL	ng/mL Range
Control A	1:50	20	17 – 23
Control B	1:50	10	8,2 – 11,8
Control C	1:50	4.1	3.4 – 4.8

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