HER-2/neu ELISA

REF 06489876

WILEX

IVD

INSTRUCTIONS FOR USE

MODE D'EMPLOI

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Cambridge, MA 02140 USA

Intended Use

The HER-2/neu ELISA is intended to quantitatively measure HER-2/neu protein in serum of women with metastatic breast cancer. The use of the HER-2/neu ELISA is indicated for follow-up and monitoring of patients with metastatic breast cancer whose initial serum HER-2/neu value is greater than 15 ng/mL. HER-2/neu values should be used in conjunction with information available from clinical and other diagnostic procedures in the management of metastatic breast cancer. The clinical utility of the serum measurement of HER-2/neu as a prognostic indicator for early recurrence and in the management of patients on immunotherapy regimens has not been fully established.

FOR IN VITRO DIAGNOSTIC USE

Method Principle	Solid Phase Sandwich ELISA
Analytical Range	0 ng/mL to 35 ng/mL
Specimen Type	Human Serum
Sample Test Volume	10.0 microliters
Sensitivity	1.5 ng p105/mL serum
Purchase of this kit licenses its use under the following patents: US patents 5,401,638 and 5,604,107; EPO patents 0494135 and 0412116; Canadian patent 2,026,250-8	

Table of Contents

Intended Use	2	Assay Procedure	4
Background	2	Sample Standard Curve	4
Principle of the Assay	2	Evaluation of Results	4
Summary of Procedure	2	Assay Quality Control	5
Standards Traceability	2	Expected Values	5
Materials Provided	2	Nonclinical Performance Characteristics	5
Safety and Warning Phrases	3	Troubleshooting	6
Materials Required but Not Provided	3	Stability and Storage Information	6
Precautions and Recommendations	3	Technical Support	6
Sample Preparation	3	References	51
Detailed Protocol	3	Explanation of Symbols	53

Background

The HER-2/neu oncogene, also referred to as c-erbB-2, encodes a protein with a molecular weight of 185,000 Daltons (p185) and belongs to a family of epithelial growth factor receptors structurally related to the human epidermal growth factor receptor [1]. The full-length p185 HER-2/neu protein is composed of a cytoplasmic domain with tyrosine kinase activity, a transmembrane domain, and an extracellular domain (ECD) that is shed from the surface of breast cancer cells [2,3]. Numerous studies have shown that the shed ECD of HER-2/neu is a glycoprotein with a molecular weight between 97 and 115 kDa and is designated p105 [3,4]. The ECD can be quantified accurately in serum with an assay [4] that uses monoclonal antibodies [5] that are directed to the external epitopes of the HER-2/neu protein. Many publications show that the ECD is shed into the blood of normal individuals and can be measured at high levels in women with metastatic breast cancer [4,6–26]. Many of these serum HER-2/neu studies have confirmed the substantial data from tissue studies that increased expression of HER-2/neu is a marker of poor prognosis, shorter overall survival, and biological aggressiveness.

Recent scientific studies suggest that quantitation of the ECD may have several important clinical applications such as monitoring breast cancer patients with metastatic disease and monitoring breast cancer patients for early recurrence [6–10,12,15,17–20,22–26]. These reports have shown that 30–50% of women with positive HER-2/neu tumors at primary diagnosis develop elevated levels of serum HER-2/neu with progression to metastatic breast cancer [4,6–26]. These studies also have illustrated that monitoring serum ECD levels post-surgery correlated with clinical course of disease and that serum HER-2/neu levels were observed to increase with disease progression or to decrease with response to therapy [12,15,19,20,25,26]. Several reports also show that elevated levels of serum HER-2/neu can occur in women with metastatic breast cancer who had primary breast tumors that were HER-2/neu negative by immunohistochemistry [6,11,18,21,27]. According to many immunohistochemistry and serum studies, the HER-2/neu protein is overexpressed in many tumors of epithelial origin including lung [28], prostate [29], pancreatic [30], colon [31], stomach [32], ovarian [33], and hepatocellular cancers [34].

Patents have been granted related to the quantitation and detection of the ECD p105 domain (in the US, patent #5,401,638; in Canada #2,026,250-8; and in Europe #0494135) as well as patents related to the quantitation and detection of the full-length p185 molecule (in the US, patent #5,604,107 and in Europe, patent #0412116). Similar patents are pending in Japan.

Principle of the Assay

The HER-2/neu ELISA is a sandwich enzyme immunoassay that utilizes a mouse monoclonal antibody for capture and a different biotinylated mouse monoclonal antibody for the detection of human HER-2/neu protein. Both capture and detector reagents specifically bind to the extracellular domain of HER-2/neu protein. The Capture Antibody has been immobilized on the interior surface of microtiter plate wells. To perform the test, an appropriate volume of specimen is incubated in the coated well to allow binding of the antigen by the Capture Antibody. The immobilized antigen is then reacted with the detector antiserum. The amount of Detector Antibody bound to antigen is measured by binding it with a streptavidin/horseradish peroxidase Conjugate, which then catalyzes the conversion of the chromogenic Substrate o-phenylenediamine (OPD) into a colored product. The colored reaction product is quantitated by spectrophotometry and is related to the amount of HER-2/neu protein in the sample.

For instructions, see the Detailed Protocol and Evaluation of Results sections of this booklet.

Summary of Procedure

Steps	Incubations
1. Add samples and Standards to wells	3 hours, 37°C
2. Wash	
3. Add Detector Antibody to wells	1 hour, 37°C
4. Wash	
5. Add Conjugate Antibody to wells	30 minutes, RT*
6. Wash	
7. Add Substrate to wells	45 minutes, RT*
8. Add Stop Solution to wells	
9. Read plate at 490 nm	

*Room temperature

Standards Traceability

The HER-2/neu ELISA was developed and calibrated by WILEX Inc., Cambridge MA, 02140, USA.

The HER-2/neu ELISA standards are calibrated against a master calibrator stock material maintained in Cambridge MA. This material is identified as HER-2/neu Master Calibrator lot R5345, and is a series of six (6) liquid reagents, spanning the range of zero to 35 ng/mL HER-2/neu p105 protein. The mass determinations for the HER-2/neu Master Calibrator have been assigned based on assay of the calibrator material in direct comparison to three (3) highly purified HER-2/neu p105 protein samples. Mass equivalents (moles) of the purified samples were determined by quantitative amino acid analysis.

Materials Provided

The following components are supplied.

Microtiter plate—Precoated microtiter plate supplied ready to use, with 96 wells (12 strips of eight) in foil zip-lock bags with a desiccant pack. Wells are coated with monoclonal anti-HER-2/neu protein antibody.

HER-2/neu Standards—Six (6) separate vials of recombinant HER-2/neu p105. Standards are calibrated in ng/mL and are labeled with values that are 50-fold greater than the actual vial dosage. Assigning these label values to a standard curve obviates the need to correct the reported dose for a 1:50 diluted sample (2% serum in buffer). See Evaluation of Results for details.

Standard#	ng/mL	Volume
6	35.0	1 mL
5	25.0	1 mL
4	15.0	1 mL
3	7.5	1 mL
2	2.5	1 mL
1	0.0	1 mL

Sample Diluent—One (1) bottle containing BSA and 0.09% sodium azide.

Detector Antibody—One (1) bottle supplied ready to use, containing biotinylated mouse monoclonal anti-HER-2/neu protein antibody in 0.01 M PBS (pH 7.4), protein stabilizer, and 0.09% sodium azide.

Conjugate Diluent—One (1) bottle containing 0.01 M PBS (pH 7.4), BSA, and 0.1% chloroacetamide.

Conjugate Concentrate—One (1) vial containing 50X streptavidin/horseradish peroxidase in buffer. Must be diluted to 1X with Conjugate Diluent to make Working Conjugate. See Table 1.

Substrate Diluent—0.1 M citrate buffer (pH 5.0) and 0.01% H₂O₂ (hydrogen peroxide).

Substrate—One (1) vial containing OPD tablets. Must be dissolved in Substrate Diluent (1 tablet/4 mL) to make Working Substrate. See Table 1.

Stop Solution—One (1) bottle containing 2.5 N H₂SO₄ (sulfuric acid).

Plat wash Concentrate (20X)—One (1) bottle. Dilute one (1) part concentrate in 19 parts high-quality water prior to use.

Safety and Warning Phrases

Harmful!

R22 – Harmful if swallowed.

S28 – After contact with skin, wash immediately with plenty of water.

Contains Sodium Azide

(Standards, Sample Diluent, Detector Antibody)

Harmful!

R40 – Possible risks of irreversible effects.

R43 – May cause sensitization by skin contact.

S36/37 – Wear suitable protective clothing and gloves.

Contains ortho-Phenylenediamine dihydrochloride

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Irritant!

R36/38 – Irritating to eyes and skin.

S36/37/39 – Wear suitable protective clothing, gloves, and eye/face protection.

S26 – In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Contains Sulfuric Acid

(Stop Solution)

Irritant!

R43 – May cause sensitization by skin contact.

S24 – Avoid contact with skin.

S37 – Wear suitable gloves.

Contains Chloroacetamide

(Conjugate Diluent and Conjugate Concentrate)

Materials Required but Not Provided

- Dry heat incubator capable of maintaining 37°C
- Pipettors: 2–20 µL, 20–200 µL, and 200–1000 µL pipettors with disposable tips
- Precision repeating pipettor
- Automated 96-well microtiter plate washer
- 12 x 75 mm culture tubes for sample preparation
- Vortex mixer
- Microtiter plate reader capable of measuring absorbance in 96-well plates at a wavelength of 490 nm
- 500- or 1000-mL graduated cylinder
- Reagent reservoirs
- Deionized water
- Plastic wrap or adhesive plate sealers
- Liquid household bleach for inactivating clinical specimens and decontamination of plate washer
- Disposable paper towels

Controls—HER-2/neu ELISA Controls consisting of recombinant p105 in Sample Diluent are available from WILEX Inc. Refer to HER-2/neu ELISA Controls, Part No. 06489884, when ordering. Controls are diluted at 1:50 prior to assay. Volumes are 0.5 mL each. Store between 2–8°C.

Precautions and Recommendations

- Store components at 2–8°C. Do not expose reagents to excessive light. Do not freeze any of the kit components.
- Do not use kit reagents past the expiration date.
- Use only the microtiter wells provided with the kit.
- Rinse all detergent residue from glassware.
- Use deionized water of the highest quality.
- Do not mix reagents from different kits.
- The buffers and reagents used in this kit contain either sodium azide or chloroacetamide as preservatives. Care should be taken to avoid direct contact with these reagents.
- Do not mouth pipet or ingest any of the reagents.
- Do not smoke, eat, or drink when performing the assay or where samples or reagents are handled.
- Human samples may be contaminated with infectious agents. Do not ingest, expose to open wounds, or breathe aerosols. Wear protective gloves and dispose of biological samples properly.
- Do not handle the Substrate Tablets with fingers or permit contact with skin, metal, or oxidizing agents. Dispose of OPD-containing solutions in compliance with local regulations.
- Wear disposable gloves and eye protection when handling Stop Solution (2.5 N sulfuric acid).
- Dispose of all working solutions (Plat wash, Conjugate, Substrate) at the end of each day.
- Prepare fresh working solutions for each subsequent assay.

CAUTION: This device contains material of human or animal origin and should be handled as a potential carrier and transmitter of disease.

Sample Preparation

For samples (and Controls), an initial dilution of 1:50 into Sample Diluent must be prepared prior to assay. Remove any flocculent material from samples by microcentrifugation prior to dilution. Label tubes and add 0.98 mL of Sample Diluent to each, followed by 0.02 mL of specimen/Control. Mix gently and avoid foaming of the diluted sample. Diluted samples and Controls may be stored in capped tubes for up to one week at 2–8°C, and six months at –20°C. Protect from light exposure. Re-assay of some samples at a higher dilution may be required. Discard any 2–8°C-stored samples showing signs of contamination. Redilute from serum.

Detailed Protocol

RECOMMENDED PROCEDURES

1. Addition of reagents must be in the order specified.
2. All six (6) Standards and test specimens should be run in duplicate.
3. Equilibrate all reagents to room temperature (15–30°C) before use.
4. Dispense only enough of each assay reagent for the number of wells being used. For every 8-well strip, dispense 1 mL of Detector Antibody, Conjugate, etc.
5. Prepare a plate map to use as a guide for location of samples, Standards, and Control wells.
6. Transfer of samples, Standards, and Controls from dilution tubes into wells can be simplified by using semi-automated, multi-channel pipettors. Contact Technical Service for recommended devices.
7. Preparation of Plat wash
 - a. If the Plat wash Concentrate is cold, allow it to reach room temperature (15–30°C) before use. Make sure all crystals are dissolved. If desired, warm at 37°C and mix well.
 - b. Dilute one (1) volume of Plat wash Concentrate with 19 volumes of distilled or deionized water. Mix well. This solution is Plat wash. The total volume required will depend on the washing method/instrument used. Approximately 1 L of this solution is required to prime an automated washer and run one microtiter plate; about 700 mL is required for each microtiter plate when manual washing is performed. Plat wash must be freshly prepared each day. Do not store Plat wash.
 - c. Plat wash must be freshly prepared each day. Do not store Plat wash.
8. When using manual methods, use caution when inverting the microtiter plate to decant or blot; press the side tabs of the frame inward to prevent the strips from falling out.
9. Microtiter plate washing is best performed by an automated plate washer. Plate washing equipment must be properly adjusted, cleaned, and maintained. Automated washers with 96 ports are recommended. Dummy 8-well strips can be used to fill the unused portion of the holder for assays using a partial plate. Retain and label previously used strips for this purpose.

- Set the fill volume to 300 μL /well. Prime the instrument with Platemash. For 96-port washers, use two 3-cycle washes. After the initial 3-cycle wash, rotate the plate 180 degrees and repeat. If a strip washer is used, perform a single wash cycle across the plate, and repeat five more times.
- After the final wash, invert the microtiter plate and strike it on an absorbent surface. Visually check that all wells are empty.


Assay Procedure

- Remove the microtiter plate from the bag. From the number of specimens to be tested and replicate determinations desired (duplicates recommended), calculate the number of strips needed. (Each Standard requires two wells; one well is needed for Substrate blank determination.) Store unused strips in the zip-lock bag with desiccant at 2–8°C (see Stability and Storage Information).
- Dilute specimens and Controls 1:50 (2%) with Sample Diluent.
- Mix the Standards, Controls, and specimens thoroughly and add 100 μL to duplicate wells. Set up one well with 100 μL of Sample Diluent for use as Substrate blank. Set up duplicate wells for each of the Controls.
- Cover the plate (food wrap or adhesive plastic cover). Incubate for 3 hours at 37°C.
- Carefully remove the plastic cover and wash the microtiter plate with Platemash.
- Add 100 μL of Detector Antibody to all wells except the Substrate blank well. Cover and incubate at 37°C for 1 hour.
- During the incubation with Detector Antibody, prepare Working Conjugate by diluting the Conjugate Concentrate with Conjugate Diluent. See Table 1 for the quantities to use for the number of strips being run.
- Wash the microtiter plate with Platemash.
- Add 100 μL of Working Conjugate to all wells except the Substrate blank well. Cover and incubate at room temperature (15–30°C) for 30 minutes.
- During the incubation with Working Conjugate, prepare Working Substrate by dissolving Substrate Tablets in Substrate Diluent. See Table 1 for the quantities to use for the number of strips being run. Vortex vigorously to assure complete dissolution. Once prepared, Working Substrate should be used within 30 minutes. Avoid exposure to light.
- Wash the microtiter plate with Platemash.
- Including the Substrate blank well, add 100 μL of Working Substrate to all wells. Cover with a fresh plastic cover and incubate the plate in the dark at room temperature (15–30°C) for 45 minutes.
- Add 100 μL of Stop Solution to each well to stop the reaction.
- Read the absorbance at 490 nm within 30 minutes.

TABLE 1. HER-2/neu ASSAY—PREPARATION OF ASSAY REAGENTS

# Strips Used	Conj. Concentrate	Conj. Diluent	Substrate Tablets	Substrate Diluent
1	20 μL	0.98 mL	1	4 mL
2	40 μL	1.96 mL	1	4 mL
3	60 μL	2.94 mL	1	4 mL
4	80 μL	3.92 mL	1	4 mL
5	100 μL	4.90 mL	2	8 mL
6	120 μL	5.88 mL	2	8 mL
7	140 μL	6.86 mL	2	8 mL
8	160 μL	7.84 mL	2	8 mL
9	180 μL	8.82 mL	3	12 mL
10	200 μL	9.80 mL	3	12 mL
11	220 μL	10.78 mL	3	12 mL
12	240 μL	11.76 mL	3	12 mL

Figure 1. Sample Standard Curve

Evaluation of Results

CONCENTRATION OF STANDARDS

The antibodies used in this assay recognize the extracellular, ligand-binding domain of the HER-2/neu protein [5]. This form of HER-2/neu has been identified with a molecular weight at or around 105 kDa. The Standards in the kit are calibrated in nanograms, which take into account this molecular weight found in serum and are prepared from a recombinant form of the 105 kDa fragment of HER-2/neu.

CONCENTRATION OF UNKNOWNNS

- Average the absorbance values for each Standard, Control, and specimen dilution to obtain the mean absorbances.
- On graph paper, plot the mean absorbance for each Standard on the y-axis versus the concentration of HER-2/neu protein (ng/mL) on the x-axis and connect points.
- Determine the concentration of HER-2/neu protein for each specimen dilution by interpolation from the standard curve. Software packages are available (such as SoftMax Pro™, Molecular Devices, Sunnyvale, CA; KC4™, Bio-Tek Instruments, Inc., Winooski, VT) which can simplify this process. Use a quadratic curve fit (second order polynomial) algorithm.
- Note: Do not assign “blank” wells using software. This will subtract the average blank readings from all other wells. It is useful for quality control and troubleshooting purposes to be able to inspect the absorbance values reported for all wells without adjustments applied to the raw data.
- Results for samples are expressed in nanograms per mL by reading directly from the standard curve values (ng/mL) as designated on the vials and in the Materials Supplied section of this booklet. For convenience, no mathematical dilution correction is needed for 1:50 diluted samples since the actual concentration in the Standard preparations are at 2% of the labeled dosage (i.e., they have been prediluted at 1:50).

DILUTION OF SAMPLES WITH HIGH CONCENTRATIONS

- For samples that give absorbance (OD) values exceeding the range of the standard curve, subsequent assay at greater dilutions will be necessary.
- To prepare further dilutions, always begin with an initial 1:50 dilution (see Sample Preparation), and then dilute serially 1:2 into Sample Diluent. Example:

Sample Dilution	Volume of Previous Dilution	Volume of Sample Diluent
1:100	0.5 mL of 1:50	0.5 mL
1:200	0.5 mL of 1:100	0.5 mL
1:400	0.5 mL of 1:200	0.5 mL

- If further diluted samples produce OD values of < 0.3, they must be re-assayed using less dilute material. The Dilution Correction Factor (see Step 4 below), when multiplied by results from such a low OD value, may result in an erroneously high HER-2/neu protein estimation.
- Any results from further diluted samples will require correcting the values obtained from the assay for any dilution beyond 1:50. Example:

Sample Dilution	Dilution Correction Factor (multiply reported result by)
1:100	2
1:200	4
1:400	8

Assay Quality Control


It is recommended that each assay meet the following performance parameters:

- Examine the results for the Controls run in the assay to confirm that the recoveries fall within the expected ranges (noted in the product insert) for each level.
- Inspect the dynamic range of the standard curve. Level 1 standard should fall between 0.04 to 0.09 OD and Level 6 standard should fall between 1.5 to 2.4 OD.
- If curve fitting software is used, inspect the calculation for R². This should fall between 0.997 and 1.0. If any of these parameters are not met, consideration should be made to repeat sample testing in a subsequent assay.

Expected Values HEALTHY INDIVIDUALS

As with all tests, each laboratory should establish its own reference range. In a population of 241 healthy women, 95% of the serum HER-2/neu protein values were found to be less than 13.7 ng/mL. There was no significant difference in HER-2/neu values between pre- and post-menopausal women. Based on this same population, the upper limit of normal (mean + 2SD) was 14.7 ng/mL. ROC analysis confirms 15 ng/mL as the appropriate cut-off between normal and elevated serum HER-2/neu protein levels. The distribution of serum HER-2/neu levels for the entire population (n = 241) is shown in Figure 2.

Figure 2. Distribution of serum HER-2/neu protein in a population of healthy women.

MONITORING OF METASTATIC BREAST CANCER PATIENTS

The clinical utility of the HER-2/neu ELISA in longitudinal monitoring of metastatic breast cancer patients was evaluated using retrospective serum samples from patients with late-stage breast cancer, covering a 6- to 12-month period, at three clinical sites in the United States. Fifty-six patients undergoing therapy whose initial serum HER-2/neu oncoprotein level was elevated (15 ng/mL or greater) were evaluated for correspondence of changes in their serum HER-2/neu levels with changes in their clinical course of disease.

A visit-to-visit analysis of the study results is presented in Table 2. The changes in serum HER-2/neu level from one visit to the next were calculated for each patient. These changes were separated into two groups: Group I had changes in serum HER-2/neu which paralleled the clinical course of disease, and Group II had changes in serum HER-2/neu level that did not parallel the clinical course of disease. Clinical course or status was determined by the physician. Correspondence of HER-2/neu changes to clinical status was determined as follows: An increase of 20% or greater from the previous visit was reflective of disease progression. If the change was less than a 20% increase from the previous visit, this was reflective of a lack of disease progression during therapy (including responding or stable disease). Stable and responding were consolidated since both reflect effectiveness of a current therapy. The 20% criterion was derived from the longitudinal variability in normal patients (determined from the serum HER-2/neu levels in six serial samples from each of 38 healthy women).

TABLE 2. CORRESPONDENCE OF METASTATIC BREAST CANCER PATIENT DISEASE STATUS WITH CHANGES IN SERUM HER-2/neu LEVELS: VISIT-BY-VISIT

Change in HER-2/neu	Change in Clinical Status		Total
	Progression	No Progression	
≥ 20% increase	52	33	85
< 20% increase	55	96	151
Total	107	129	236

Concordance = 62.7% (CI = 56.4 to 68.6%)

Predictive Value: Lack of Progression = 63.6% (CI = 55.7 to 70.8%)

Progression = 61.2% (CI = 50.6 to 70.8%)

The following figures show typical examples of changes in serum HER-2/neu levels over time for two of the 56 metastatic breast cancer patients from the clinical study. Figure 3 shows a patient whose disease is responding and Figure 4 shows a patient with progressive disease.


Figure 3. Monitoring of a 38-year-old Stage IV breast cancer patient with the HER-2/neu ELISA. Longitudinal changes in serum HER-2/neu levels correlate with changes in disease status.


Figure 4. Monitoring of a 74-year-old Stage IV breast cancer patient with the HER-2/neu ELISA. Longitudinal changes in serum HER-2/neu levels correlate with changes in disease status.

Nonclinical Performance Characteristics

IMPRECISION

Imprecision was determined by testing four control samples consisting of three buffer-based Controls and one human serum pool. Three different kit lots were used at each of three different laboratories. Triplicate determinations per run for each Control were determined in 130 runs. Within-run and total precision were calculated for each Control based on the NCCLS EP-5A protocol. Within-run imprecision for three sites and kit lots were at or below 10% CV with total imprecision across all kits and sites between 10 and 17.7% CV.

Control Sample	# of Runs	# of Replicates	Mean (ng/mL)
Control 1	130	385	24.5
Control 2	131	390	9.8
Control 3	131	386	3.3
Serum 4	130	385	9.5

Control Sample	Within-Run Precision		Total Precision	
	Std Dev	% CV	Std Dev	% CV
Control 1	1.46	6.0	2.62	10.7
Control 2	0.69	7.0	1.03	10.4
Control 3	0.34	10.2	0.59	17.7
Serum 4	0.65	6.9	1.02	10.8

SENSITIVITY

Minimal detectable concentration of analyte was determined by repeated measurement of a zero dose sample (Sample Diluent). Mean and standard deviation of apparent recovery of HER-2/neu were calculated for 230 determinations made over a period of 20 days between three different kit lots and three different laboratories. After addition of two standard deviations to the mean recovered dose of HER-2/neu, the minimal detectable concentration was determined to be 1.5 ng of p105/mL serum.

ANALYTICAL SPECIFICITY

Cross-reactivity. Epidermal growth factor receptor (EGFr) is a closely related molecule to the HER-2/neu growth factor receptor at 88% homology with the intracellular domain and 44% homology with the extracellular domain.

We challenged the HER-2/neu ELISA with EGFR at 850 ng/mL. This material produced no signal in the assay, indicating that the immunoreagents in the device produce no cross-reaction with EGFR.

Interfering substances. Human anti-mouse antibody (HAMA), found in rare instances in some human serum, has the potential to bind to the critical reagents in immunoassays, causing false positive or negative signals. Several known HAMA-positive serum samples as well as rheumatoid-factor-positive samples were tested in the device. No false signals were generated with these samples. Device reagents have been formulated to block any such false reactivity.

Serum HER-2/neu measurements might be performed while patients are taking vitamins, over-the-counter drugs, or undergoing chemotherapy. To test the possibility that such agents could interfere with accurate HER-2/neu determinations, potential exogenous interferents were spiked into a Positive Control serum containing a known concentration of HER-2/neu and were subsequently tested in the assay. Some common serum constituents also were tested for their potential to act as endogenous interferents. None of the compounds tested had an effect on analyte recovery. See Table 3.

ANALYTICAL RANGE

The HER-2/neu ELISA is capable of accurate quantitation of serum HER-2/neu levels in prediluted samples, which provide an assay response in the range of 1.5 to 35 ng of p105 HER-2/neu protein per mL serum. Samples that produce a signal in excess of the upper limit of the standard curve range (35 ng/mL) must be further diluted with Sample Diluent and re-tested in the assay. Be certain to correct the reported recovery of HER-2/neu protein for dilution preparations that are greater than the routine 1:50 dilution.

TABLE 3. POTENTIAL ELISA INTERFERENTS

POTENTIAL INTERFERENTS	TEST CONCENTRATIONS	POTENTIAL INTERFERENTS	TEST CONCENTRATIONS
ENDOGENOUS		OTC DRUGS, ETC.	
Triglycerides (Intralipin-20%)	3000.0 mg/dL	Acetaminophen	200.0 µg/mL
Hemoglobin	1.0 g/dL	Aspirin	500.0 µg/mL
Immunoglobulin (gamma globulin)	6.0 g/dL	Ibuprofen	400.0 µg/mL
Albumin	6.5 g/dL	Caffeine	100.0 µg/mL
Bilirubin	25.0 mg/dL		
Heparin	0.46 mg/mL		
EXOGENOUS		CHEMOTHERAPEUTIC AGENTS	
Vitamin A trans Retinol Acetate	10.0 IU/mL	Aminoglutethamide	398.0 µg/mL
Vitamin B1 Thiamine	3.0 µg/mL	Bleomycin	0.16 U/mL
Vitamin B2 Riboflavin	3.4 µg/mL	Cis-Platin	173.0 µg/mL
Vitamin B6 Pyroxide HCL	4.0 µg/mL	Doxorubicin	51.8 µg/mL
Vitamin B12	12.0 ng/mL	Cyclophosphamide	800.0 µg/mL
Vitamin C Ascorbic Acid	30.0 µg/mL	Diethylstilbestrol	23.0 µg/mL
Vitamin D2 Ergocalciferol	0.8 IU/mL	Estramustine	102.2 µg/mL
Vitamin E Tocopherol Acetate	0.06 IU/mL	Flutamide	10.0 µg/mL
Folic Acid	0.8 µg/mL	5-Fluorouracil	1600.0 µg/mL
Niacin	40.0 µg/mL	Lupron	15.0 µg/mL
		Methotrexate	450.0 µg/mL
		Mitoxantrone	56.0 µg/mL
		Mitomycin C	73.0 µg/mL
		Megace	27.0 ng/mL
		Tamoxifen	1.4 µg/mL
		Vincristine	4.88 µg/mL
		Vinblastine	16.3 µg/mL
		Herceptin (Trastuzumab)	500 µg/mL

LIMITATIONS OF THE PROCEDURE

In performing steps in manual immunoassays, there is the potential of erroneous measurement (sample recovery) as a result of variation in the incubation time that samples are subjected to in each of the incubation steps. All wells should receive essentially the same incubation time required for each reagent step. In particular, when adding Standards, Controls, and samples to wells in the plate, the time between addition of the first sample and the last should be within, at most, a 15-minute time frame. Prepare all samples in advance. Use of semi-automated, multi-tip pipettors is highly recommended for shortening the times for these proce-

dures. In order to prevent cross-contamination, be sure to use fresh tips for each sample addition.

Immunoassays using so-called sandwich configurations, such as this device, are subject to interference and false results arising from samples that contain immunoreactive agents against mouse antibodies. HAMA (human anti-mouse antibody), heterophilic antibody, and rheumatoid factor are the primary examples of these interferents. Precautions have been taken by the manufacturer in formulating reagents contained in this device to minimize or eliminate such interference. However, care should be taken when evaluating assay results that may be inconsistent with the overall clinical status of the patient in light of this possible interference.

It is possible that a patient with confirmed breast cancer may have serum HER-2/neu protein levels within the range of those observed in healthy individuals. Use caution in the interpretation of HER-2/neu results. Serum HER-2/neu levels may not always indicate the presence or absence of malignant disease. The HER-2/neu value should be used as part of an overall clinical assessment that includes additional clinical evaluation and diagnostic tests.

Troubleshooting

- Each assay must include the six (6) Standards tested in duplicate using the protocol described in the Detailed Protocol. Incubation times or temperatures significantly different from those specified may give erroneous results.
- The shape of the standard curve is nonlinear and varies slightly from assay to assay. For data reduction by hand, apply point-to-point graphing. If using ELISA software, apply a quadratic curve fit (second order polynomial) algorithm.
- Poor duplicates, if accompanied by elevated values for the zero Standard, indicate insufficient washing. If all instructions in the Detailed Protocol were followed accurately, such results indicate a need for washer maintenance. Also, avoid introducing bubbles into wells when pipetting samples and reagents.
- The Substrate blank well should read less than or equal to 0.05 absorbance units. A possible cause of higher values is exposure of the Working Substrate to light either before or during the incubation step. Working Substrate should be used within 30 minutes of preparation.
- Low signal may indicate 1) that drying of the plate has occurred between steps. Do not allow the plate to dry out. Add the next reagent immediately after wash; or 2) that improper washing has been performed. Be sure the microtiter plate washer has been properly primed with Platetwash (See Detailed Protocol).

Stability and Storage Information

All of the reagents included in the HER-2/neu ELISA have been tested for stability. Reagents should not be used beyond the stated expiration date. Kit reagents should be stored at 2–8°C, with the exception of Platetwash Concentrate, which may be stored at room temperature. Assay plates should be stored sealed in the original foil bag with desiccant pack. Once opened, the device will perform within specifications for up to 4 weeks (30 days) within the expiration date period indicated on the label.

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References

1. Coussens L, et al. Tyrosine kinase receptor with extensive homology to EGF receptor shares chromosomal location with neu oncogene. *Science* 1985;230:1132–1139.
2. Brandt-Rauf P, Pincus MR, Carney WP. The c-erbB-2 protein in oncogenesis: Molecular structure to molecular epidemiology. *Critical Reviews in Oncogenesis* 1994;5(2&3):313–329.
3. Zabrecky JR, et al. The extracellular domain of p185/neu is released from the surface of human breast carcinoma cells, SK-BR-3. *Journal of Biological Chemistry* 1991;266(3):1716–1720.
4. Carney WP, et al. Detection and quantitation of the human neu oncoprotein. *Journal of Tumor Marker Oncology* 1996;(2):53–72.
5. McKenzie SJ, et al. Generation and characterization of monoclonal antibodies specific for the human neu oncogene product, p185. *Oncogene* 1989;4:543–548.
6. Andersen TI, et al. Detection of c-erbB-2 related protein in sera from breast cancer patients. Relationship to erbB-2 gene amplification and c-erbB-2 protein. *Acta Oncologica* 1995;34(4):499–504.
7. Watanabe N, et al. Serum c-erbB-2 in breast cancer patients. *Acta Oncologica* 1994;33(8):901–904.
8. Molina R, et al. c-erbB-2 oncoprotein in the sera and tissue of patients with breast cancer. Utility in prognosis. *Anticancer Research* 1996;16(4B):2295–2300.
9. Molina R, et al. Utility of c-erbB-2 in tissue and in serum in the early diagnosis of recurrence in breast cancer patients: Comparison with carcinoembryonic antigen and CA 15.3. *British Journal of Cancer* 1996;74:1126–1131.
10. Molina R, et al. c-erbB-2 oncoprotein, CEA, and CA 15-3 in patients with breast cancer. *Breast Cancer Research & Treatment* 1998;51:109–119.
11. Fehm T, et al. The prognostic significance of c-erbB-2 serum protein in metastatic breast cancer. *Oncology* 1998;55:33–38.
12. Isola JJ, et al. Elevated erbB-2 oncoprotein levels in preoperative and follow-up serum samples define an aggressive disease course in patients with breast cancer. *Cancer* 1994;73(3):652–658.
13. Leitzel K, et al. Elevated serum c-erbB-2 antigen levels and decreased response to hormone therapy of breast cancer. *Journal of Clinical Oncology* 1995;13(5):1129–1135.
14. Hayes DF, et al. Elevated circulating c-neu oncogene product in patients with breast cancer. Abstract presented at the 12th Annual San Antonio Breast Cancer Symposium, 1989. (December 8–9, 1989.)
15. Kath R, et al. The neu-oncogene product in serum and tissue of patients with breast carcinoma. *Annals of Oncology* 1993;4:585–590.
16. Yamauchi H, et al. Prediction of response to anti-estrogen therapy in advanced breast cancer patients by pre-treatment circulating levels of extracellular domain of the HER-2/neu protein. *Journal of Clinical Oncology* 1997;15(7):2518–2525.
17. Mercer DW, et al. HER-2/neu oncoprotein monitoring in sera of patients with breast cancer. *Clinical Chemistry*, Abstract #846, 1995;41(S6, part 2):S226.
18. Kandl H, et al. Soluble c-erbB-2 fragment in serum correlates with disease stage and predicts shortened survival in patients with early-stage and advanced breast cancer. *British Journal of Cancer* 1994;70(4):739–742.
19. Narita T, et al. c-erbB-2 protein in the sera of breast cancer patients. *Breast Cancer Research and Treatment* 1992;24:97–102.
20. Dia J, et al. Clinical utility of an automated serum HER-2/neu assay in monitoring patients with metastatic breast cancer. *Amer Assoc Clin Chem* 1998. Abstract #208.
21. Krainer M, et al. Tissue expression and serum levels of HER-2/neu in patients with breast cancer. *Oncology*, 1997;54:475–481.
22. Kynast B, et al. Determination of a fragment of the serum c-erbB-2 translational protein p185 in serum of breast cancer patients. *J. Cancer Res Clin Oncol* 1993;119:249–252.
23. Takaaski H, et al. Clinical significance of serum c-erbB-2 protein in patients with primary breast cancer. *Oncology Reports* 1997;4(2):349–352.
24. Klein B, et al. Soluble c-erbB-2 (p185) in breast cancer patients in relation to prognosis. *Oncology Reports* 1995;2(5):759–761.
25. Imoto S, et al. Serum c-erbB-2 levels in monitoring of operable breast cancer patients. *Jpn J Clin Oncol* 1999;29(7):336–339.
26. Willsher PC, et al. Prognostic significance of serum c-erbB-2 protein in breast cancer patients. *Breast Cancer Research & Treatment* 1996;40:251–255.
27. Titus K. Another ingredient added to HER-2 mix. *CAP Today*, Nov 1999.
28. Brandt-Rauf P, et al. Detection of increased amounts of the extracellular domain of the c-erbB-2 oncoprotein in serum during pulmonary carcinogenesis in humans. *Int Journal Cancer* 1994;56:383–386.
29. Myers RB, et al. Elevated serum levels of p105 erbB-2 in patients with advanced stage prostatic adenocarcinoma. *Int Journal Cancer* 1996;69(5):398–402.
30. Okada N, et al. Elevated serum c-erbB-2 protein levels in patients with pancreatic cancer: Correlation to metastasis and shorter survival. *Oncology* 1995;52:392–396.
31. Wu JT, Astill ME, Zhang P. Detection of the extracellular domain of c-erbB-2 oncoprotein in sera from patients with various carcinomas: Correlation with tumor markers. *Journal of Clinical Laboratory Analysis* 1993;7:31–40.
32. Chariyalertsak S, et al. Comparison of c-erbB-2 oncoprotein expression in tissue and serum of patients with stomach cancer. *Tumor Biol* 1994;15:294–303.
33. Meden H, et al. Prognostic significance of p105 (c-erbB-2, HER-2/neu) serum levels in patients with ovarian cancer. *Anti Cancer Research* 1997;17:757–760.
34. Luo J-C, et al. Serum c-erbB-2 oncoprotein in hepatocellular carcinogenesis. *Medical Science Research* 1993;21:305–307.

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En: Microtiter Plate: 96 well plate pre-coated with monoclonal antibody. Serves as solid-phase capture surface for assay reaction steps.

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DILSPE

En: Sample Diluent: Reagent used to dilute specimens prior to assay.

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STD 1-6

En: Standard 1-6: Six standards contain HER-2/neu in a dose range from 0 to 35 ng/mL in neutral buffer. Standards are used to calibrate each assay and determine levels of HER-2/neu in samples.

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please contact WILEX Inc.