

Siemens Medical Solutions Diagnostics 511 Benedict Avenue Tarrytown, NY 10591; USA	Emergency Telephone 1-800-424-9300 (Chemtrec) Technical Information 1-877-229-3711
--	---

1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND THE COMPANY/UNDERTAKING

Product Name: TIMP-1 ELISA

Product Number: 06489973

Date Prepared: 3/27/2007

Revision Number: 2

2. COMPOSITION/INFORMATION ON INGREDIENTS

This product contains several components. Please refer to the individual Safety Data Sheets for the components listed below.

<i>Component Name</i>	<i>Product Number</i>
TIMP-1 Microtiter Plate	26-7010
Sample Diluent	26-7020
TIMP-1 Standards, Levels 1-6	26-7025-1/6
Conjugate	26-7045
Platewash Concentrate	26-7095
Stop Solution	26-7092
Filtered RO Water	26-7000
Substrate (Liquid pNPP)	26-7053

Products may contain components which require reporting under SARA Title III Section 313 (40 CFR 372). See applicable component MSDS.

EU Classification (1999/45/EC):



EU Risk and Safety Phrases: R22, S28

14. TRANSPORT INFORMATION

Proper Shipping Name: Not regulated

Technical Name:

UN Number: N/A

Hazard Class and Packaging Group: N/A

Label(s): N/A





Packing Instruction (Passenger Aircraft): N/A

Packing Instruction (Cargo Aircraft): N/A

Shipping Container Type: N/A

Sales Unit: 1 kit.

Prepared by: Regulatory Affairs Department

HMIS	Hazard	Rating	
 HEALTH		1	4 Severe Hazard
 FLAMMABILITY		0	3 Serious Hazard
 REACTIVITY		0	2 Moderate Hazard
 PERSONAL_PROTECTION		0	1 Slight Hazard
		B	0 Minimal Hazard
			* An asterisk indicates significant chronic health effects

The opinions expressed herein are those of qualified experts within Siemens Medical Solutions Diagnostics. We believe that the information contained herein is current as of the date of this Material Safety Data Sheet. Since the use of this information and these opinions and the conditions of use of the product are not within the control of Siemens Medical Solutions Diagnostics, it is the users' obligation to assure safe use of the product.